

Think Smart Inc  
SmartUpload

---

TRANSDYNE

# Installation Guide

37637 Five Mile Rd# 315, • Livonia, MI • 48154  
10225 Barnes canyon rd, Suite a208 • san diego, ca • 92121  
Phone: 866 535 3054 • Fax: 866 453 9025

SMART UPLOAD

# Installation and configuration

---

© Transdyne Inc.  
37637 Five Mile Rd • Suite 315  
Phone 866 535 3054 • Fax 866 453 9025

---

# Table of Contents

Chapter 1	-	Introduction	2
		Installation Sequence	2
		Estimated Time	2
Chapter 2	-	Requirements	3
Chapter 3	-	Begin Installation Sequence	4
Chapter 4	-	Install Certificate	10
Chapter 5	-	Starting Smart upload	12



## Introduction

*TRANSDYNE Inc., an industry leader in information systems for the healthcare industry, specializes in medical transcription and technical consulting. This guide walks you through the process of installing the software provided by Transdyne to electronically send your dictations to our servers over the internet.*

This windows based Smart upload program from TransDyne is used to download dictations from your handheld digital recorder to your local computer and then transmit them to TransDyne's servers over https using 128 bit security.

## Installation Sequence

- ü Install SMART UPLOAD
- ü Configure SMART UPLOAD

## Estimated Time

All of the above process would approximately take 10 minutes of your time, at the end of which you will be able to successfully transmit your dictations to TransDyne's servers.

## Requirements

### Hardware

Pentium II 233 MHz or Higher

64 MB Ram or Higher

1 GB Hard disk space or Higher

CD Rom Drive

### Software

OS: Windows 2000/XP/2003

Internet Explorer 6.0 or Higher.

Microsoft Word or Word Perfect for Editing.

### Administrator

You are required to login as administrator or as a user that has administrative privileges to proceed with installation. You are not required to be an administrator to use this program.

### User ID and Password


This is supplied to you in your setup kit.

### Internet Connection

## Beginning Installation process

*Smart Upload installation should start automatically as soon as you insert the CD supplied with the Install packet. If it does not you can manually start the installation by double clicking on setup.exe that is found in the install CD.*

The Install shield wizard will walk you through the installation process. Please click next to proceed with the installation on the Welcome Screen.




Customer Information Dialog: Please enter your name and organization name and select 'Anyone who uses this computer (All Users)' and click 'Next' to proceed with installation.



The image shows a Windows-style dialog box titled "Smart Upload - InstallShield Wizard". The dialog has a blue title bar with a close button (X) in the top right corner. Below the title bar, the text "Customer Information" is displayed in bold, followed by the instruction "Please enter your information." There are two text input fields: "User Name:" and "Organization:". Below these fields, the text "Install this application for:" is followed by two radio button options: "Anyone who uses this computer (all users)" (which is selected) and "Only for me ()". At the bottom left, the text "InstallShield" is visible. At the bottom right, there are three buttons: "< Back", "Next >", and "Cancel".

Destination Folder Dialog: The installation program will default the installation folder to 'C:\Program Files\TransDyne\Smart Upload'. If you wish to install to a different folder please click on the change button and select the folder that you want this program to be installed to.

Click on 'Next' to proceed with Smart Upload installation



Feature Selection Dialog: This will be defaulted to 'Complete' installation, which will install both the download from hand held feature and upload to TransDyne server feature. This program will only support some handheld recorders, please check with TransDyne's technical support team to find out if the recorder that you intend to use is supported by this program. If you intend to use this only to upload dictations to TransDyne's servers, please choose custom.

Click Next to proceed to the next step.


If you choose custom please skip the next page and go to page 9.




Please skip this step if you choose complete installation.

Smart DSS Download Installation: By default, Smart DSS Downloader will not be installed if you select custom setup.

Please click 'Next' to proceed to the next step



Ready to Install: Please click Install to begin Installation.




Click On "OK" button. You will be prompted with the following message.

## Install Certificate


This will install the certificate required by the https protocol. Please click 'Next' to proceed.





## Start Smart Upload

Check the "Launch Smart Upload" check box and click on the Finish button to start Smart Upload.




Input your user id as supplied by TransDyne in the "User Name" field.

Input your password as supplied by TransDyne in the "password" field.

Input your account id as supplied by TransDyne in the "Account ID" field.

Click 'Ok' when you are finished.



You will see the following screen when finished.



## Copyright Notice and Usage Terms

### Copyright Notice

The TransDyne SmartUpload User's Manual is Copyright © 2004- 2005 ThinkSmart Inc, U.S.A. All Rights Reserved.

### Terms of Usage for TransDyne SamrtUpload User's Manual

This manual is free for the exclusive use of TransDyne SmartUpload users. It is created to help you use the Web Delivery System. Sharing this manual with any one who is not a user of TransDyne Web Delivery is strictly prohibited.